

OLIGO NAME	REF#	CHEMISTRY	BASES	
Anchored Oligo dT (20)	UN-PR015-005	DNA	20	
SEQUENCE (5' → 3')				
5' TTT-TTT-TTT-TTT-TTT-TTT-TV 3'				
BACKBONE (BASES + LINKAGES)				
A	0 GC (%) ¹	3.0	Quantity (OD@260nm)	0.82
C	0 MW (g.mol ⁻¹) ¹	6028.4	(nmol)	5.00
G	0 Ext. coeff. (l.mol ⁻¹ .cm ⁻¹) ¹	164600	(µg)	30.1
T	19 Tm (°C) ¹	42.0		
Wobbles	1			
¹ Theoretical values				
PURIFICATION(S)				
HPLC-RP				
QUALITY CONTROL(S)				
MS - CGE				Passed

DELIVERY CONDITION
Room Temperature
RECONSTITUTION
To make up to 100 µM concentration, add 50 µl of recommended buffer. To make up to 20 µM concentration, add 250 µl of recommended buffer. To make up to 5 µM concentration, add 1000 µl of recommended buffer.
LONG TERM STORAGE
-20 °C

COMMENTS		
General: Tm (°C) calculated with 50 mM NaCl and 0 % formamide (PCR conditions) Specific:		
TECHNICAL SUPPORT		
If you have any questions feel free to call our Oligo Centre		
<table border="0" style="width: 100%;"> <tr> <td style="width: 50%; vertical-align: top;"> EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB </td> <td style="width: 50%; vertical-align: top;"> NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com </td> </tr> </table>	EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB	NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com
EUROPE ☎ 00 800 666 00 123 (European toll free number), ✉ oligo@eurogentec.com Kaneka Eurogentec S.A. Liège Science Park Rue Bois Saint-Jean 5 - 4102 SERAING BELGIUM Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88 E-mail: info@eurogentec.com Web: www.eurogentec.com RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050 BIC: BBRUBEBB	NORTH AMERICA ☎ +1 (510) 896-1350 ✉ oligo.na@egt-biotech.com EGT North America, Inc. 34801 Campus Drive - Fremont, CA 94555 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com	