

EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00
 • Fax: +32 4 372 75 00 • Toll-free: + 800 666 00 123 • info@eurogentec.com
 • www.eurogentec.com

**NORTH AMERICA**

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560
 • Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com
 • www.anaspec.com

Technical Data Sheet

iD MOPS Running Buffer Powder ID-WCRUB1-005

Eurogentec products are sold for research or laboratory use only and are not to be administered to humans or used for medical diagnostics.

Product

Cat #	Description	Size
ID-WCRUB1-005	iD MOPS Running Buffer Powder	5 packs for 5 x 1L running buffer

Description

iD MOPS Running Buffer Powder is used for iD PAGE Gels electrophoresis. Using proprietary techniques, iD MOPS Running Buffer Powder is optimized to have a long shelf life, high resolution and fast electrophoresis.

Storage

Store at RT

Reconstitution

Reconstitute with 1000 mL H₂O to make 1X running buffer per bag of powder.

For further information please contact our Customer Help Desk:**For Europe:**

E-mail: info@eurogentec.com
 Tel: +32 4 372 76 65

For USA:

E-mail: service@anaspec.com
 Tel: +1 510 791 9560