

## EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • BELGIUM • Tel.: +32 4 372 74 00  
• Fax: +32 4 372 75 00 • Toll-free: + 800 666 00 123 • info@eurogentec.com  
• www.eurogentec.com



## NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA • Tel.: +1-510-791-9560  
• Toll-free: +1 800-452-5530 • Fax: +1 510-791-9572 • service@anaspec.com  
• www.anaspec.com

## Technical Data Sheet

### iD Blot Pads ID-BLPAB1-020

Eurogentec products are sold for research or laboratory use only and are not to be administered to humans or used for medical diagnostics.

## Product

Cat #	Description	Size
ID-BLPAB1-020	iD Blot Pads (Basic)	20 packs

## Description

The iD Blot Pads are the consumable part of the iD Blot Device (discontinued) for fast electroblotting of proteins from polyacrylamide gels to nitrocellulose or PVDF membranes.

## Kit content

Component	Piece	Specifications
iD Blot Pad	20	Blotting filter paper presoaked with proprietary anode or cathode buffer
iD Blot Cathode Pad		90 mm (l) × 80 mm (w) × 2.5 mm (t)
iD Blot Anode Pad		90 mm (l) × 80 mm (w) × 2.5 mm (t)
iD Blot Equilibration Buffer (125 ml)	2	125 ml × 2
Gel Window (66 mm×56 mm)	1	Polycarbonate
Gel Window (76 mm×64 mm)	1	Polycarbonate
Gel Window (88 mm×78 mm)	1	Polycarbonate
Absorbent Filter Paper	1	80 mm × 70 mm - Vegetable fiber
Sponge Cushion	2	

## Storage

The iD Blot Pads are stable for 8 months at room temperature

## For further information please contact our Customer Help Desk:

E-mail: [info@eurogentec.com](mailto:info@eurogentec.com)

Tel: +32 4 372 76 65