

OLIGO NAME	CATALOG #	CHEMISTRY	BASES	DATE	
ROX-T10	UN-CT035-005	DNA	10	20-Jul-18	
SEQUENCE (5' → 3')					
5' ROX™ - TTT-TTT-TTT-T 3'					
BACKBONE (BASES + LINKAGES)					
A	0	GC (%) ¹	0.0	Quantity (OD@260nm)	0.52
C	0	MW (g.mol ⁻¹) ¹	3675.8	(nmol)	5.0
G	0	Ext. coeff. (l.mol ⁻¹ .cm ⁻¹) ¹	104200	(µg)	18.4
T	10	Tm (°C) ¹	20.0		
Others	1				
¹ Theoretical values					
PURIFICATION(S)					
HPLC-RP					
QUALITY CONTROL(S)					
MS					Passed

DELIVERY CONDITION

Room Temperature

RECONSTITUTION

To make up to 100 µM concentration, add 50 µl of recommended buffer.

To make up to 20 µM concentration, add 250 µl of recommended buffer.

To make up to 5 µM concentration, add 1000 µl of recommended buffer.

LONG TERM STORAGE

Store at -20 °C

COMMENTS

General: Tm % GC calculated with 50 mM NaCl and 0 % formamide (PCR conditions)

Specific: Stable until 20-Jan-20 under dry condition at room temperature, QC available on request

TECHNICAL SUPPORT

If you have any questions feel free to call our Oligo Centre

EUROPE

☎ 00 800 666 00 123 (European toll free number),

 ✉ oligo@eurogentec.com

 Kaneka Eurogentec S.A. Liège Science Park
 Rue Bois Saint-Jean 5 - 4102 SERAINING BELGIUM
 Tel: +32(0)4 372 74 00 - Fax: +32(0)4 264 07 88

 E-mail: info@eurogentec.com Web: www.eurogentec.com

 RPM Liège T.V.A.-(BE)-0427.348.346 - ING Belgique Bank - IBAN: BE86 3400 2118 6050
 BIC: BBRUBEBB

NORTH AMERICA

☎ +1 (510) 896-1350

 ✉ oligo.na@egt-biotech.com

 EGT North America, Inc.
 34801 Campus Drive - Fremont, CA 94555
 Tel.: +1 (510) 896 -1350 * Fax: +1 (510) 791-9572
 E-mail: ivd.na@egt-biotech.com * Web: www.egt-biotech.com